



Baylor College of Medicine Microarray Core Facility

One Baylor Plaza, Room 452A, Houston, Texas 77030
(713) 798-7699 microarray@bcm.tmc.edu



Silanization Protocol

1. Clean glass trays by rinsing with ddH₂O before using. The silane tray should then be rinsed with 95% EtOH, and the acetone tray should be rinsed with acetone.
2. Sonicate the slides in acetone for 10 minutes.
3. After the acetone sonication, wash slides in ddH₂O tray at least twice.
4. Sonicate the slides in .1M NaOH for 10 minutes.
5. While the slides are sonicating in NaOH, make the following silane solution under a hood:

REAGENTS	For small trays (30 slides) add:	For large trays (60 slides) add:
95% EtOH	350 mL	800 mL
Glacial Acetic Acid	42 μ L	96 μ L
Glycidoxypropyltrimethoxy silane	11 mL	25 mL

NOTE: Do not add the silane until the NaOH sonication is finished. After the silane is added, stir for 1-2 minutes.

6. After the NaOH sonication, wash slides in ddH₂O tray at least twice.
7. Sonicate the slides in silane solution for at least 3 minutes. DO NOT oversonicate (5 minutes maximum!!).
8. Wash the slides in 100% EtOH tray.
9. Dry with prepurified N₂ gas. Store in 100°C oven. Wait at least 2 hours before using slides.

2X Phosphate Printing Buffer (pH 10.5)

- 400 mM Na₂HPO₄
- 1.6 M NaCl
- NaOH to pH 10.5
- Filter Sterilized 1-2000

Add 5 μ L to printing plates and 5 μ L of cDNA to make a 1X working buffer.